



NOLAN TECHNOLOGIES LLC

---

# urTalker Pro Application



## User Guide

---

# Table of Contents

- urTalker Pro Administration ..... 1
  - Access Settings ..... 2
  - PIN ..... 2
  - Modify Categories ..... 3
  - Modify Words ..... 3
  - Modify Existing Category ..... 4
  - Add New Category ..... 4
  - Delete / Move Words ..... 5
  - Category Properties ..... 5
    - Color ..... 5
    - Sound..... 5
    - Record Sound ..... 6
  - Modify Existing Word ..... 7
  - Add New Word..... 7
- urTalker Pro User..... 10
  - Words / Images..... 11
  - Sentence Builder..... 12
- About urTalker ..... 13

---

# urTalker Pro Administration

## Initial Load Screen


When you first load the application you will get a one-time notification:

### Important Notice:

This application has an administrative feature that requires a password.

The default password is 0000. You will need this password to access the settings menu.

Press “Ok” to proceed. Then click the **settings** button in the top right to change your password or begin to setup urTalker.




## Settings – Administration - Login



### Access Settings


Click the settings button in the top right to configure urTalker Pro. This area is for adding or changing categories, words and default settings.



### Enter PIN

If you have not changed your PIN from the default setting enter 0000 otherwise enter your unique PIN in this screen

## Settings – Administration – Main Screen Options



### Loop Around

If you want image and categories to be viewable in a continuous loop, leave this box checked. If you only want images and categories to scroll one direction, uncheck this box.

### Modify Categories

To add new categories or modify existing categories select this option.

### Modify View of Categories

To set the number of categories to see at one time select this option.

### Modify Words

To add new words (images) or modify existing words select this option.


### Change PIN

To change your PIN code, select this option.

### About


Click this to find out what version of urTalker you have and to find links to our website.

# Settings – Administration – Modify Categories



## Modify Existing Category

Swipe or use the arrows at the top of the screen to find the category you want to edit. Then click the edit button located at the bottom left of the screen.




## Add New Category

To add a new category click the plus button located at the bottom left of the screen.



# Settings – Administration – Modify Categories



## Category Properties

For existing or new categories this screen will appear.

### Name

Click and type name of category. This will appear on the user screens.

### Color

To change the color select “customize color” and pick one of four defaults **or** click define “define your color” to create a color of your choosing.

### Customize Image

Click to upload or replace the category image. Note: It is not required if you just want a color.

### Sound

Sound is optional for categories. If you wish to have a sound played when a category is selected click “associate sound”.

### Default Category

To make this category the default or initial category loaded, check this box.

### Default by Schedule

To make this category appear or disappear on a day/time schedule check this box and then set the times.

### Grid

This controls the number of images to display inside this category. Click to change.

## Delete / Move Words

To delete a category click the “delete” button. Caution: if you delete a category and the words have not been reassigned those words will not be moved to another category.

If any words are assigned, please use the “move all words to another category” button before deleting. This will safely reassign those words to another category.

## Settings – Administration – Sound



### Record Sound



Click this button to record a sound using your device.


### Use iTunes Library

If you already have a sound click this option and select it from your library.





### Record Audio

To record new audio click the  button to start recording and then click the  button to stop recording.



### Review and Save

To review your recorded audio click . To delete and re-record click .

Once you are ready to save your sound click “associate sound”




# Settings – Administration – Modify Words



## Modify Existing Word

Swipe or use the arrows at the top of the screen to find the word you want to edit. Then click the edit button located at the bottom left of the screen.



## Add New Word

To add a new word click the plus button located at the bottom left of the screen.




# Settings – Administration – Add New Words



## Step 1: Name & Priority

Enter the name for the word and then select a priority for the category. Ex. If you want this word to appear first select 1. If you don't select they will sort alphabetical. If 2 words have the same priority then they will sort alphabetical.



## Step 2: Take Photo / Select Image

Click the take photo or library option



## Select Image from Roll


Select the image you want to use



## Crop/Positioning

You can zoom in or out as well as reposition the image if needed.

## Settings – Administration – Add New Words (cont'd)




### Step 3: Associate Sound

Select the option for adding sound.




### Record Sound

See steps in modify category for recording sound. Once complete click “associate sound” button.



### Step 4: Assign Category

Click the appropriate option of what category the word should go into.



### Save Word

Click “save word” to save and complete the add word process.

## urTalker Pro User

### Initial Screen - Categories



In user mode the initial launch will start with a set number of categories as seen. This configuration, modify view of categories, is set at 9. (Note: you can have more than nine and the categories will then swipe or click to view more. The starting category is based on your default setting.

Swipe left or right to move through the categories or use the buttons located at the top of the screen.

Click the category to view words/images inside of that category.

# urTalker Pro User : Word / Image Playback


## Words / Images

Once inside a category all the words/images available will be displayed based on the category view setting. This category is set to nine.


Swipe left or right to view additional words or use the arrows at the bottom of the screen.

Click the middle button to go back to categories or click the home button at the top right to go back to the initial category screen.

Press any image to play the sound for that image/word.




# urTalker User : Sentence Builder Playback



## Sentence Builder


To build a sentence, click the “+” sign in the bottom right and it will add that word to your sentence builder.

Continue pressing the “+” sign and add additional words to your sentence.

To play the entire sentence click the play  button located in the bottom right.

To remove a word from your sentence click the “-“ button.

To re-arrange your sentence just click and drag the images.



## About urTalker

urTalker is all about enabling individuals with a voice.

Our applications are designed with the individual and their unique needs in mind. While every application comes with default categories and images, every solution provides the ability to create any number or additional categories, load any type of image and most importantly record your voice for use with those images and words.

Our Purpose- urTalker was developed by the founders to help their non-verbal autistic child, have a voice. Many non-verbal, autistic and special needs individuals have the ability to learn and use assistive technology, but the price for specific devices can be prohibitive. urTalker was developed to be affordable .

Our Future- this is just the beginning. We are open to feedback and will continue to update and improve our apps for your benefit or the benefit of your loved one. Please feel free to contact us or contribute to the product so we continue to provide the best communication solution possible.